

ACT tubes available

The MAX-ACT™ and the C-ACT™ activated clotting tubes are now available from Gribbles Veterinary*.

- Provides an easy assessment of Activated Clotting Time (ACT)
- Results in minutes
- Easy to operate
- Only one reagent
- Sample volumes of 0.5mL (MAX-ACT™) or 2mL (C-ACT™)
- Room temperature storage

C- ACT™ tubes requires 2 mls blood
Cost: \$3.25 + GST /tube

MAX-ACT™ tubes requires 0.5mL blood
Cost: \$4.95 + GST/tube

Minimum order: five tubes per pack

For further details contact your local CSM or
our Help Desk on 1300 307 190.

Reference: *AM See, KL Swindwells, MJ Sharman, KL Haack, D Goodman, A Delaporta, I Robertson and SF Foster. Activated coagulation times in normal cats and dogs using MAX-ACT™ tubes.
Aust. Vet Journal, Vol 87, No. 7, July 2009.

Fax: 03 9538 6741
E-mail: vets@gribbles.com.au



Gribbles
VETERINARY

